Babi.Plus® Bubble PAP^{MB} Valve

Accurate system for delivery of BCPAP therapy

The Babi.Plus Bubble PAP^{MB} Valve provides a convenient method of delivering continuous PAP (positive airway pressure) to infants weighing less than 10 kg. The patented design of the Bubble PAP^{MB} Valve and the Multi-Bubble drop tube combine macro and micro bubble oscillations to deliver accuracy and stability throughout the infant's course of therapy, freeing the clinician to focus on patient care.

The Babi.Plus Bubble PAP^{MB} Valve was the first to receive FDA 510k market clearance for delivery of positive end expiratory pressure in patients weighing < 10 kg in hospital environments to increase end lung pressure above atmospheric in a constant flow condition. The PAP is adjustable from 1 to 10 cm H_2O, using gas flows from 1 to 12 liters per minute for precise pressure delivery. The Bubble PAP^{MB} Valve's ergonomic design allows pressure to be set easily and prevents unintended changes in PEEP levels.

The Babi.Plus Bubble PAP^{MB} Valve jar has an integrated ISO compliant equipment mount, which makes it easy to secure the jar onto a standard equipment pole with a clamp. The Venti.Plus™ Single or Dual Pole Clamp companion products are designed to pair with Babi.Plus products and work well with other products requiring pole clamps.

The Babi.Plus Bubble PAP^{MB} Valve offers a simpler, less invasive therapy that supports gas exchange across large and small airways through the application of the MultiBubble. It can be used in all neonatal critical care units, delivery rooms, pediatric critical care units and special procedure units.

Ordering Information

1002	Babi.Plus BCPAP Kit	12/case
1000	Babi.Plus Bubble PAP^{MB} Valve 0-10 cm H_2O *	12/case

Companion Products

1010	Babi.Plus nCPAP Kit, Small
1020	Babi.Plus nCPAP Kit, Large
1030	Babi.Plus Pressure Relief Manifold
B020	Babi.Plus nCPAP Nasal Cannula, Sz 0
B021	Babi.Plus nCPAP Nasal Cannula, Sz 1
B022	Babi.Plus nCPAP Nasal Cannula, Sz 2
B023	Babi.Plus nCPAP Nasal Cannula, Sz 3
B024	Babi.Plus nCPAP Nasal Cannula, Sz 4
B025	Babi.Plus nCPAP Nasal Cannula, Sz 5
B026	Babi.Plus nCPAP Nasal Cannula, Sz 6
B027	Babi.Plus nCPAP Nasal Cannula, Sz 7

1005	Venti.Plus Single Pole Clamp
1006	Venti.Plus Dual Pole Clamp
30100	Baby Nose Bumper™ and RespiraGel Mustache™
30150	Baby Nose Bumper Replacement
30300 Series	Baby Cap™ and Circuit Bumpers™ in seven sizes
30400	Baby Chin Strap™
30500 Series	Baby Head Band™ and Circuit Bumpers in four sizes
30510	Baby Band™ for Chin and Head

* US patent, with additional patents pending



FEATURES AND BENEFITS

- Accurate method to deliver PAP for infants weighing ≤10 kg
- Pressure settings clearly marked for precise pressure adjustment
- Positive engagement pressure lever reduces risk of inadvertent changes
- 360° swivel inlet connector accepts 15 mm OD by 22 mm ID expiratory circuit connector
- Water Level adjustment port allows fluid to be added or removed without patient disconnect
- Single-patient use closed system minimizes the incidence of atmospheric fluid contamination
- Latex free
- Ready for immediate set up in all clinical environments

Respiralogics products are available from Respiralogics and our specialty respiratory care distributors. For our latest product information visit www.respiralogics.com or email us at 4info@respiralogics.com.

Respiralogics is a Global Respiratory Solutions, Inc. Company
3545 Airway Drive, Suite 104 • Reno, NV 89511 • 775 954 0160 • www.respiralogics.com

© 2026 Respiralogics. All rights reserved. Respiralogics and the Respiralogics logo, Baby Cap and Circuit Bumpers, Baby Nose Bumper and RespiraGel Mustache, Baby Chin Strap, Baby Head Band and Circuit Bumpers, and Baby Bands for Chin and Head are trademarks of Global Respiratory Solutions, Inc. Babi.Plus is a registered trademark and Venti.Plus is a trademark of Galemed Corporation. C1000.7