

NeoTree™ Tube Support

Unique holder adapts to neonatal
and small pediatric circuits

respiralogics™

NeoTree Tube Support works in conjunction with neonatal/small pediatric ventilator circuits. The Tube Support is reversible in position to secure tubes for neonatal/small pediatric patients.


NeoTree assists in the reduction of tension on the endotracheal tube or face mask to prevent kinking of the circuit, inadvertent extubations and/or movement of the endotracheal tube or mask.

NeoTree is made of lightweight aluminum with a baked-on powder coat finish. The tube holder is durable and easily cleaned. NeoTree Tube Support provides a stable holder for tubes in all positions.

Ordering Information for NeoTree Tube Support

420-1948 NeoTree Tube Support (1/box)

Respiralogics products are available from Respiralogics and our specialty respiratory care distributors. For our latest product information visit www.respiralogics.com or email us at information@respiralogics.com.



FEATURES

- Reversible holder to secure term, pre-term neonatal and small pediatric breathing circuits
- Constructed of lightweight aluminum
- Multiple cut-outs for ease of placement
- Easy to clean powder coating finish
- Reusable and durable
- May be used in NICU, PICU, surgery and clinics for stable tube positioning

Respiralogics is a Global Respiratory Solutions, Inc. Company
3545 Airway Drive, Suite 104 • Reno, NV 89511 • 775 954 0160 • www.respiralogics.com

© 2017 Respiralogics. All rights reserved. Respiralogics and the Respiralogics logo are trademarks of Global Respiratory Solutions, Inc. NeoTree is a trademark used under license by McArthur Medical Sales Inc. C1948.1