



ODAPT Adapter for soft mist inhalers (SMI) is the only product created specifically for use with the soft mist family of inhaled medications. Designed by a respiratory therapist to provide clinicians and patients a solution for face mask and tracheostomy application of soft mist inhalers.

ODAPT aids in the delivery of SMI medication for pediatric, adult and elderly patients where a face mask is required. It also creates a new option for delivery of SMI medication for all patients with a tracheostomy.

ODAPT Adapters are latex, BPA, and DEHP free and interface with commercially available face masks. Available as adapter only or with face mask. Don't adapt it, ODAPT it!



### Ordering Information for ODAPT adapters

420-7190	ODAPT adapter for soft mist inhaler (40/box)
420-7193	ODAPT adapter for soft mist inhaler with adult mask, yellow (10/box)
420-7194	ODAPT adapter for soft mist inhaler with large adult mask, green (10/box)

Respiralogics products are available from Respiralogics and our specialty respiratory care distributors. For our latest product information visit [www.respiralogics.com](http://www.respiralogics.com) or email us at [information@respiralogics.com](mailto:information@respiralogics.com).



Respiralogics is a Global Respiratory Solutions, Inc. Company  
3545 Airway Drive, Suite 104 • Reno, NV 89511 • 775 954 0160 • [www.respiralogics.com](http://www.respiralogics.com)

© 2017 Respiralogics. All rights reserved. Respiralogics and the Respiralogics logo are trademarks of Global Respiratory Solutions, Inc. ODAPT is a trademark used under license by McArthur Medical Sales Inc, and holds registered Canadian and international design patents. Respimat is a registered trademark of Boehringer Ingelheim International GmbH, used under license. C7190.1



### ODAPT ADAPTERS

Cost-effective solution to allow RESPIMAT® soft mist inhalers to be used with standard masks

Allows patient to self-administer inhaled medications via mask or tracheotomy tube

Standard fitting allows direct connection to commercially available face mask or tracheostomy tube

### ODAPT WITH MASK

Easier to coordinate administration of soft mist inhalers

Ideal for patients who have difficulty creating a strong seal around SMI mouthpiece

Anatomically shaped mask with soft seal lip for optimum seal

Allows consistency for delivery of SMI medication by mask (patient preference)