

# Babi.Plus® nMask

Silicone resuscitation mask for optimal comfort and seal

respiralogics™

Babi.Plus nMask is precisely designed to comfortably conform to an infant's face, facilitating a good seal for optimal ventilation and non-invasive support.



The circular design of the nMask allows easy creation of a comfortable anatomical seal for the very small premature infant to toddler face. The clear nMask provides a good view of the patient/mask interface to ensure proper fit and size.

The silicone nMask is available in six sizes; ELBW through small pediatric patient. Each nMask is individually packaged and ready for use with resuscitation bag or T-piece resuscitator system.

The Babi.Plus nMask can be used in all neonatal critical care units, delivery rooms, pediatric critical care units and special procedure units.

## Ordering Information

ANJ0001	nMask 25 mm (1 piece/bag, 5 bags/box)
ANJ0002	nMask 30 mm (1 piece/bag, 5 bags/box)
ANJ0003	nMask 35 mm (1 piece/bag, 5 bags/box)
ANJ0004	nMask 40 mm (1 piece/bag, 5 bags/box)
ANJ0005	nMask 45 mm (1 piece/bag, 5 bags/box)
ANJ0006	nMask 50 mm (1 piece/bag, 5 bags/box)

Respiralogics products are available from Respiralogics and our specialty respiratory care distributors. For our latest product information visit [www.respiralogics.com](http://www.respiralogics.com) or email us at [information@respiralogics.com](mailto:information@respiralogics.com).



**FEATURES**

- Designed specifically to conform to premature and newborn faces
- Ergonomic design with soft material to provide better seal and comfort
- Silicone material to minimize skin irritation; bio-compatible, hypoallergenic and latex free
- Available in six sizes; fits ELBW infant through small pediatric faces
- 15mm connector with two-level barb design for easier attachment to resuscitator devices or T-piece resuscitator circuit
- Single patient use
- Individually packaged and easy to apply to neck of resuscitation instrument

Respiralogics is a Global Respiratory Solutions, Inc. Company  
3545 Airway Drive, Suite 104 • Reno, NV 89511 • 775 954 0160 • [www.respiralogics.com](http://www.respiralogics.com)

© 2017 Respiralogics. All rights reserved. Respiralogics and the Respiralogics logo are trademarks of Global Respiratory Solutions, Inc. Babi.Plus is a registered trademark of Galemed Corporation. ANJ0001.1