

# Babi.Plus® Pressure Relief Manifold

Prevent pressure build up in the circuit

respiralogics™

12.5 cm H<sub>2</sub>O Pressure Relief Manifold with Gas Supply Tube and Oxygen Sensor Port is a single use device for use with patients weighing ≤ 10 kg. The manifold is designed to prevent breathing circuit pressure from exceeding the preset pressure limit.

The Babi.Plus Pressure Relief Manifold also includes a Gas Supply Tube, an Oxygen Sensor Port and capped port for an oxygen sensor analyzer.

The Babi.Plus Pressure Relief Manifold is designed to pair with the Babi.Plus Bubble PAP Valve 0-10 cm H<sub>2</sub>O and Babi.Plus Cannula Kits companion products. The Babi.Plus Pressure Relief Manifold can be used in neonatal critical care units, delivery rooms, pediatric critical care units and special procedure units.

## Ordering Information

1030	Babi.Plus Pressure Relief Manifold (12/case)
1002	Babi.Plus BCPAP Kit (12/case)

## Companion Products

1010	Babi.Plus Bubble PAP Valve	30100	Baby Nose Bumper and RespiraGel Mustache
1010	Babi.Plus nCPAP Kit, Small	30301	Baby Cap & Circuit Bumpers, Size 1
1020	Babi.Plus nCPAP Kit, Large	30302	Baby Cap & Circuit Bumpers, Size 2
1016	Babi.Plus nCPAP Complete Kit, Small	30303	Baby Cap & Circuit Bumpers, Size 3
1026	Babi.Plus nCPAP Complete Kit, Large	30304	Baby Cap & Circuit Bumpers, Size 4
1030	Babi.Plus Pressure Relief Manifold	30305	Baby Cap & Circuit Bumpers, Size 5
B020	Babi.Plus nCPAP Nasal Cannula, Sz 0	30306	Baby Cap & Circuit Bumpers, Size 6
B021	Babi.Plus nCPAP Nasal Cannula, Sz 1	30307	Baby Cap & Circuit Bumpers, Size 7
B022	Babi.Plus nCPAP Nasal Cannula, Sz 2	30400	Baby Chin Strap
B023	Babi.Plus nCPAP Nasal Cannula, Sz 3		
B024	Babi.Plus nCPAP Nasal Cannula, Sz 4		
B025	Babi.Plus nCPAP Nasal Cannula, Sz 5		
B026	Babi.Plus nCPAP Nasal Cannula, Sz 6		
B027	Babi.Plus nCPAP Nasal Cannula, Sz 7		
1005	Venti.Plus Single Pole Clamp		
1006	Venti.Plus Dual Pole Clamp		

Respiralogics products are available from Respiralogics and our specialty respiratory care distributors. For our latest product information visit [www.respiralogics.com](http://www.respiralogics.com) or email us at [information@respiralogics.com](mailto:information@respiralogics.com).



## FEATURES

Easy to use manifold for delivery of non-invasive therapy for infants weighing ≤ 10 kg

15 mm ID x 22 mm OD inlet adaptor interfaces with all humidifier chambers

Integrated oxygen supply tubing with fits all connector for mixed gas delivery

Constructed of Phthalate Free, Latex free and hypoallergenic material

Pre-assembled, single patient use manifold

Respiralogics is a Global Respiratory Solutions, Inc. Company  
3545 Airway Drive, Suite 104 • Reno, NV 89511 • 775 954 0160 • [www.respiralogics.com](http://www.respiralogics.com)

© 2017 Respiralogics. All rights reserved. Respiralogics and the Respiralogics logo are trademarks of Global Respiratory Solutions, Inc. Babi.Plus is a registered trademark and Venti.Plus is a trademark of Galemed Corporation. C1030.4