

Babi.Plus® BCPAP Kit

respiralogics™

Bubble PAP Valve and Pressure Relief Manifold kit for easy-to-use, safe and effective BCPAP therapy

The Babi.Plus BCPAP Kit provides a convenient method for delivering continuous PAP (positive airway pressure) to infants weighing less than (<) 10 kg by combining the patented Babi.Plus Bubble PAP Valve and Babi.Plus Pressure Relief Manifold in a single easy-to-assemble package.

The Bubble PAP Valve and the Pressure Relief Manifold are conveniently packaged for delivery of PAP in patients weighing < 10 kg. The Bubble PAP Valve is adjustable from 1 to 10 cm H₂O, using gas flows from 1 to 12 liters per minute for precise pressure delivery. The Pressure Relief Manifold with gas supply tube and capped oxygen sensor port is designed to prevent breathing circuit pressure from exceeding the preset pressure limit.

The BCPAP Kit is a single-use device designed to allow ease of set up without the cumbersome and time-consuming tasks associated with bubble devices. Babi.Plus BCPAP Kit's universal connections attach to standard humidifiers, delivery circuits, Babi.Plus Cannula Kits, complementary Baby products and other commercially available nasal interfaces.

Ordering Information

1002	Babi.Plus BCPAP Kit	12/case
------	---------------------	---------

Companion Products

1010	Babi.Plus nCPAP Kit, Small	30100	Baby Nose Bumper™ and RespiraGel Mustache™
1020	Babi.Plus nCPAP Kit, Large	30150	Baby Nose Bumper Replacement
B020	Babi.Plus nCPAP Nasal Cannula, Sz 0	30301	Baby Cap™ & Circuit Bumpers™, Size 1
B021	Babi.Plus nCPAP Nasal Cannula, Sz 1	30302	Baby Cap & Circuit Bumpers, Size 2
B022	Babi.Plus nCPAP Nasal Cannula, Sz 2	30303	Baby Cap & Circuit Bumpers, Size 3
B023	Babi.Plus nCPAP Nasal Cannula, Sz 3	30304	Baby Cap & Circuit Bumpers, Size 4
B024	Babi.Plus nCPAP Nasal Cannula, Sz 4	30305	Baby Cap & Circuit Bumpers, Size 5
B025	Babi.Plus nCPAP Nasal Cannula, Sz 5	30306	Baby Cap & Circuit Bumpers, Size 6
B026	Babi.Plus nCPAP Nasal Cannula, Sz 6	30307	Baby Cap & Circuit Bumpers, Size 7
B027	Babi.Plus nCPAP Nasal Cannula, Sz 7	30400	Baby Chin Strap™
1005	Venti.Plus™ Single Pole Clamp	30501	Baby Head Band™ & Circuit Bumpers, XS
1006	Venti.Plus Dual Pole Clamp	30502	Baby Head Band & Circuit Bumpers, S
GB10	GiO™ 1 Pressure Gauge	30503	Baby Head Band & Circuit Bumpers, M
		30504	Baby Head Band & Circuit Bumpers, L
		30510	Baby Band™ for Chin and Head

Respiralogics products are available from Respiralogics and our specialty respiratory care distributors. For our latest product information visit www.respiralogics.com or email us at 4info@respiralogics.com.



FEATURES AND BENEFITS

- Complete kit ready for nCPAP delivery
- Pressure control settings clearly marked for easy and secure adjustment of nCPAP
- 360° swivel inlet connector accepts standard 15 mm OD by 22 mm ID circuit fittings
- Monitor and adjust fluid level without interrupting BCPAP delivery
- Pressure Relief Manifold connects directly to gas supply, humidifier and oxygen analyzer sensor
- Closed system minimizes spillage and atmospheric contamination of H₂O
- Constructed of phthalate-free, latex-free, hypoallergenic material

Respiralogics is a Global Respiratory Solutions, Inc. Company
3545 Airway Drive, Suite 104 • Reno, NV 89511 • 775 954 0160 • www.respiralogics.com

© 2019 Respiralogics. All rights reserved. Respiralogics, the Respiralogics logo, Baby Nose Bumper and RespiraGel Mustache, Baby Head Band and Circuit Bumpers, Baby Chin Strap, Baby Cap and Circuit Bumpers, and Baby Band for Chin and Head are trademarks of Global Respiratory Solutions, Inc. Venti.Plus and GiO are trademarks and Babi.Plus is a registered trademark of GaleMed Inc. C1002.1