Babi.Plus® BCPAP Kit

Bubble PAP^{MB} Valve and Pressure Relief Manifold kit for easy-to-use, safe, and effective BCPAP therapy

respiralogics™

The Babi.Plus BCPAP Kit provides a convenient method for delivering continuous PAP (positive airway pressure) to infants weighing less than (<) 10 kg. The kit combines the patented design of the Babi.Plus Bubble PAP^{MB} Valve* with the Babi.Plus Pressure Relief Manifold for a single easy-to-assemble package.

The Bubble PAP^{MB} Valve and the Pressure Relief Manifold are conveniently packaged for delivery of PAP in patients weighing < 10 kg. The Bubble PAP^{MB} Valve is adjustable from 1 to 10 cm H_2O, using gas flows from 1 to 12 liters per minute for precise pressure delivery. The Pressure Relief Manifold with gas supply tube and capped oxygen sensor port is designed to prevent breathing circuit pressure from exceeding the preset pressure limit.

The BCPAP Kit is a single-use device designed to allow ease of set up without the cumbersome and time-consuming tasks associated with bubble devices. Babi.Plus BCPAP Kit's universal connections attach to standard humidifiers, delivery circuits, Babi.Plus Cannula Kits, complementary Baby products and other commercially available nasal interfaces.

Ordering Information

1002	Babi.Plus BCPAP Kit	12/case

Companion Products

1010	Babi.Plus nCPAP Kit, Small
1020	Babi.Plus nCPAP Kit, Large
B020	Babi.Plus nCPAP Nasal Cannula, Sz 0
B021	Babi.Plus nCPAP Nasal Cannula, Sz 1
B022	Babi.Plus nCPAP Nasal Cannula, Sz 2
B023	Babi.Plus nCPAP Nasal Cannula, Sz 3
B024	Babi.Plus nCPAP Nasal Cannula, Sz 4
B025	Babi.Plus nCPAP Nasal Cannula, Sz 5
B026	Babi.Plus nCPAP Nasal Cannula, Sz 6
B027	Babi.Plus nCPAP Nasal Cannula, Sz 7

1005	Venti.Plus™ Single Pole Clamp
1006	Venti.Plus Dual Pole Clamp
30100	Baby Nose Bumper™ and RespiraGel Mustache™
30150	Baby Nose Bumper Replacement
30300 series	Baby Cap™ & Circuit Bumpers™, seven sizes
30400	Baby Chin Strap™
30500 series	Baby Head Band™ & Circuit Bumpers, four sizes
30510	Baby Band™ for Chin and Head

Respiralogics products are available from Respiralogics and our specialty respiratory care distributors. For our latest product information visit www.respiralogics.com or email us at 4info@respiralogics.com.

* US patent, with additional patents pending



FEATURES AND BENEFITS

Complete kit ready for nCPAP delivery in all clinical environments

Pressure control settings clearly marked for easy and secure adjustment of nCPAP

360° swivel inlet connector accepts standard 15 mm OD by 22 mm ID circuit fittings

Monitor and adjust fluid level without interrupting BCPAP delivery

Pressure Relief Manifold connects directly to gas supply, humidifier and oxygen analyzer sensor

Closed system minimizes spillage and atmospheric contamination of H_20

Constructed of phthalate-free, latex-free, hypoallergenic material

Respiralogics is a Global Respiratory Solutions, Inc. Company
3545 Airway Drive, Suite 104 • Reno, NV 89511 • 775 954 0160 • www.respiralogics.com

© 2026 Respiralogics. All rights reserved. Respiralogics, the Respiralogics logo, Baby Nose Bumper and RespiraGel Mustache, Baby Head Band and Circuit Bumpers, Baby Chin Strap, Baby Cap and Circuit Bumpers, and Baby Band for Chin and Head are trademarks of Global Respiratory Solutions, Inc. Venti.Plus is a trademark and Babi.Plus is a registered trademark of GaleMed Inc. C1002.3