Babi.Plus® nCPAP Kit

Safe, comfortable method to deliver nCPAP and NIV

respiralogics™

The Babi.Plus nCPAP Kit is a single use device for use with neonates, infants and small pediatric patients weighing ≤ 10 kg to deliver intermittent or continuous gas flow via a nasal prong. The contoured nasal prongs made with medical grade silicone minimize skin irritation and provide stability during non-invasive therapy. The unique cannula body reduces contact or "grabbing" of the infant's skin while keeping the cannula comfortably in place. An integrated pressure monitoring line allows monitoring of proximal airway pressure without having lines or ports near the infant's face.

Offered in eight different sizes, Babi.Plus cannulas can be used on a wide range of patients. The Small Babi.Plus nCPAP Kit is packaged with four sizes of silicone nasal prongs individually packaged for infants ≤ 2 kg. The Large Babi.Plus nCPAP Kit is packaged with four sizes of silicone prongs individually packed for infants ≥ 2 kg up to 10 kg.

The Cannula Kit contains a prong sizing guide, dual-limb circuit to connect cannula, and two silicone/rubber adapters to connect the distal ends of the dual-limb circuit to the humidifier delivery circuit.

The Babi.Plus Bubble PAP^{MB} Valve, Babi.Plus BCPAP Kit, Baby Cap™ and Circuit Bumpers™, Baby Nose Bumper™, Baby Chin Strap™, the Venti.Plus™ Single Pole Clamp and Venti.Plus Dual Pole Clamp are companion products designed to pair with the Cannula Kit. The Babi.Plus Cannula Kit can be used in neonatal critical care units, delivery rooms, pediatric critical care units and special procedure units.

Ordering Information

1010	Babi.Plus nCPAP Kit, Small: Size 0, 1, 2 and 3	12/case
1020	Babi.Plus nCPAP Kit, Large: Size 4, 5, 6 and 7	12/case

Companion Products

1002	Babi.Plus BCPAP Kit
1000	Babi.Plus Bubble PAP^{MB} Valve 0-10 cm H_2O*
1030	Babi.Plus Pressure Relief Manifold
1005	Venti.Plus Single Pole Clamp
1006	Venti.Plus Dual Pole Clamp
30100	Baby Nose Bumper and RespiraGel Mustache™
30150	Baby Nose Bumper Replacement
30400	Baby Chin Strap

* US patent, with additional patents pending

30301	Baby Cap & Circuit Bumpers, Size 1
30302	Baby Cap & Circuit Bumpers, Size 2
30303	Baby Cap & Circuit Bumpers, Size 3
30304	Baby Cap & Circuit Bumpers, Size 4
30305	Baby Cap & Circuit Bumpers, Size 5
30306	Baby Cap & Circuit Bumpers, Size 6
30307	Baby Cap & Circuit Bumpers, Size 7
30501	Baby Head Band™, Size Extra Small
30502	Baby Head Band, Size Small
30503	Baby Head Band, Size Medium
30504	Baby Head Band, Size Large
30510	Baby Band™ for Chin and Head



FEATURES AND BENEFITS

- Comfortable method for delivery of non-invasive therapy via nasal prongs for infants weighing less than or equal to 10 kg
- Nasal prongs are made of medical grade silicone to minimize skin irritation
- Contoured base of prong allows better fit without resting on philtrum
- Nasal prongs are designed to keep prong body away from septum and prevent pinching of septum
- Constructed of DEHP-free, latex-free and hypoallergenic material
- Both single-patient-use kits packaged with 4 cannulas, providing replacement cannulas in one package
- Integrated pressure line for proximal airway pressure measurement

Respiralogics products are available from Respiralogics and our specialty respiratory care distributors. For our latest product information visit www.respiralogics.com or email us at 4info@respiralogics.com.

Respiralogics is a Global Respiratory Solutions, Inc. Company
3545 Airway Drive, Suite 104 • Reno, NV 89511 • 775 954 0160 • www.respiralogics.com

© 2026 Respiralogics. All rights reserved. Respiralogics and the Respiralogics logo, Baby Cap and Circuit Bumpers, Baby Nose Bumper and RespiraGel Mustache, Baby Chin Strap, Baby Head Band and Circuit Bumpers, and Baby Bands for Chin and Head are trademarks of Global Respiratory Solutions, Inc. Babi.Plus is a registered trademark and Venti.Plus is a trademark of Galemed Corporation. C1010.12