Babi.Plus®
Pressure Relief Manifold

Prevent pressure build up in the circuit

12.5 cm H_2O Pressure Relief Manifold with Gas Supply Tube and Oxygen Sensor Port is a single use device for use with patients weighing ≤ 10 kg. The manifold is designed to prevent breathing circuit pressure from exceeding the preset pressure limit.

The Babi.Plus Pressure Relief Manifold also includes a Gas Supply Tube, an Oxygen Sensor Port and capped port for an oxygen sensor analyzer.

The Babi.Plus Pressure Relief Manifold is designed to pair with the Babi.Plus Bubble PAP^{MB} Valve 0-10 cm H_2O and Babi.Plus Cannula Kits companion products. The Babi.Plus Pressure Relief Manifold can be used in neonatal critical care units, delivery rooms, pediatric critical care units and special procedure units.

Ordering Information

1030	Babi.Plus Pressure Relief Manifold	12/case
1002	Babi.Plus BCPAP Kit	12/case

Companion Products

1000	Babi.Plus Bubble PAP^{MB} Valve
1010	Babi.Plus nCPAP Kit, Small
1020	Babi.Plus nCPAP Kit, Large
B020	Babi.Plus nCPAP Nasal Cannula, Sz 0
B021	Babi.Plus nCPAP Nasal Cannula, Sz 1
B022	Babi.Plus nCPAP Nasal Cannula, Sz 2
B023	Babi.Plus nCPAP Nasal Cannula, Sz 3
B024	Babi.Plus nCPAP Nasal Cannula, Sz 4
B025	Babi.Plus nCPAP Nasal Cannula, Sz 5
B026	Babi.Plus nCPAP Nasal Cannula, Sz 6
B027	Babi.Plus nCPAP Nasal Cannula, Sz 7
1005	Venti.Plus™ Single Pole Clamp
1006	Venti.Plus Dual Pole Clamp

30100	Baby Nose Bumper™ and RespiraGel Mustache™
30150	Baby Nose Bumper Replacement
30301	Baby Cap & Circuit Bumpers™, Size 1
30302	Baby Cap & Circuit Bumpers, Size 2
30303	Baby Cap & Circuit Bumpers, Size 3
30304	Baby Cap & Circuit Bumpers, Size 4
30305	Baby Cap & Circuit Bumpers, Size 5
30306	Baby Cap & Circuit Bumpers, Size 6
30307	Baby Cap & Circuit Bumpers, Size 7
30400	Baby Chin Strap™
30501	Baby Head Band™, Size Extra Small
30502	Baby Head Band, Size Small
30503	Baby Head Band, Size Medium
30504	Baby Head Band, Size Large
30510	Baby Band™ for Chin and Head



FEATURES AND BENEFITS

Easy to use manifold for delivery of non-invasive therapy for infants weighing ≤ 10 kg

15 mm ID x 22 mm OD inlet adaptor interfaces with all humidifier chambers

Integrated oxygen supply tubing with fits-all connector for mixed gas delivery

Constructed of Phthalate Free, Latex free and hypoallergenic material

Pre–assembled, single patient use manifold

Respiralogics products are available from Respiralogics and our specialty respiratory care distributors. For our latest product information visit www.respiralogics.com or email us at 4info@respiralogics.com.

Respiralogics is a Global Respiratory Solutions, Inc. Company
3545 Airway Drive, Suite 104 • Reno, NV 89511 • 775 954 0160 • www.respiralogics.com

© 2026 Respiralogics. All rights reserved. Respiralogics and the Respiralogics logo, Baby Cap and Circuit Bumpers, Baby Nose Bumper and RespiraGel Mustache, Baby Chin Strap, Baby Head Band and Circuit Bumpers, and Baby Bands for Chin and Head are trademarks of Global Respiratory Solutions, Inc. Babi.Plus is a registered trademark and Venti.Plus is a trademark of Galemed Corporation. C1030.8