Babi.Plus® nTestLung

Neonatal Lung Simulator

respiralogics™

The Babi.Plus nTest Lung is an easy to use neonatal test lung for demonstration and testing applications with infant and neonatal ventilator support devices. The nTest Lung consists of two individual silicone lungs and a Y connector. The nTest Lung is the ideal instrument for training and demonstration in addition to performing ventilator verification testing.

Specifications

Description	Specification
Compliance	Approximate 2.0 mL/cm H2O
Tidal Volume	25 ml each silicone bag, 50 ml total
Resistance	300 cm H2O/L/s at 1 L/sec
Interface Connector	Y connector with 15/10 mm port
Straight Connector	10 mm OD
Pressure Monitoring Port	4.3 mm ID
Weight	43 grams

Ordering Information

2010	Babi.Plus nTest Lung	1 each

Respiralogics products are available from Respiralogics and our specialty respiratory care distributors. For our latest product information visit www.respiralogics.com or email us at 4info@respiralogics.com.



FEATURES AND BENEFITS

- Two lung feature with Y connector to mimic anatomical design
- Standard ISO connections
- Lightweight, easy to use and portable
- Manufactured of high performance material
- Packaged in convenient zippered canvas bag for storage and transport

Respiralogics is a Global Respiratory Solutions, Inc. Company
3545 Airway Drive, Suite 104 • Reno, NV 89511 • 775 954 0160 • www.respiralogics.com

© 2026 Respiralogics. All rights reserved. Respiralogics and the Respiralogics logo are trademarks of Global Respiratory Solutions, Inc.
Babi.Plus is a registered trademark of Galemed Corporation. C2010.6