Baby Chin Strap™

For babies who cannot keep their mouths closed

Baby Chin Strap is a single-patient-use device intended to help keep small patients' mouths closed during delivery of nCPAP and NIV therapy. Baby Chin Strap offers support with a soft, skin-friendly strap that cups the chin and is secured to the Baby Head Band or Baby Cap with hook and loop tabs. Baby Chin Strap is adjustable to provide support to easily close the mouth and provides a comfortable solution for mouth leaks.

Baby Chin Strap may be used in conjunction with Baby Head Band™ and Circuit Bumpers™ or Baby Cap™ and Circuit Bumpers. Designed for use in neonatal critical care units, delivery rooms, pediatric critical care units and special procedure units for stable delivery of BCPAP of NIV.





FEATURES AND BENEFITS

- Intended for very low birth weight through full-term infants
- Assists in keeping the mouth closed to optimize nCPAP and NIV therapy
- Split chin design allows for ease of application and cupping the chin
- Made of skin friendly, soft knit to optimize infant comfort
- Single patient use, latex free

Ordering Information

30400	Baby Chin Strap	20 / box

Companion Products

30301	Baby Cap and Circuit Bumpers, Size 1
30302	Baby Cap and Circuit Bumpers, Size 2
30303	Baby Cap and Circuit Bumpers, Size 3
30304	Baby Cap and Circuit Bumpers, Size 4
30305	Baby Cap and Circuit Bumpers, Size 5
30306	Baby Cap and Circuit Bumpers, Size 6
30307	Baby Cap and Circuit Bumpers, Size 7

Respiralogics products are available from Respiralogics and our specialty respiratory care distributors. For our latest product information visit www.respiralogics.com or email us at 4info@respiralogics.com.

Respiralogics is a Global Respiratory Solutions, Inc. Company
3545 Airway Drive, Suite 104 • Reno, NV 89511 • 775 954 0160 • www.respiralogics.com

© 2026 Respiralogics. All rights reserved. Respiralogics, the Respiralogics logo, Baby Head Band and Circuit Bumpers, Baby Chin Strap and Baby Cap and Circuit Bumpers are trademarks of Global Respiratory Solutions, Inc. C30400.2